

SPECTRUM VASCULAR Arterial PRESSURE MONITORING SETS & TRAYS



We've Got You Covered.

With over 40 years as a recognized leader in minimally invasive technologies, you can rely on us. Our portfolio is diverse and is designed to meet the challenges you face.



www.spectrumvascular.com

WHY SPECTRUM VASCULAR

Used for arterial pressure monitoring and blood sampling.^{†,††,†††}

FEATURES & BENEFITS

- Inserted via the minimally invasive Seldinger insertion technique to minimize vessel trauma.
- Lines available for radial or femoral insertion to increase placement options.
- Femoral access set is available with extra-firm nylon catheter.

PRIMARY USE AREAS

- Intensive Care Unit
- Perioperative
- Emergency Room



~1 in 6
individuals hospitalised with
influenza had severe outcomes¹

Intended for arterial blood pressure monitoring and blood sampling in adult and pediatric patients.

[†] 2.5 French catheters are intended for patients from birth and older.

^{††} 3.0 and 4.0 French catheters are intended for patients 1 year and older.

^{†††} 5.0 French catheters are intended for patients aged 12 years and older.



¹Sumner, et al. Severity of influenza-associated hospitalisations by influenza virus type and subtype in the USA, 2010-19: a repeated cross-sectional study, *The Lancet Microbe*, Volume 4, Issue 11, 2023, Pages e903-e912, ISSN 2666-5247, [https://doi.org/10.1016/S2666-5247\(23\)00187-8](https://doi.org/10.1016/S2666-5247(23)00187-8).

www.spectrumvascular.com

Available in a range of sizes and design configurations.



Product Name	Order Number	Reference Part Number (RPN)	Catheter Fr/ Equivalent gage	Material	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Tip Configuration	Needle Baseplate
Sets								
Arterial Pressure Monitoring Set	G58400	APM-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
Arterial Pressure Monitoring Set	G58401	APM-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
Arterial Pressure Monitoring Set	G58402	APM-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
Arterial Pressure Monitoring Set	G58403	APM-2.5-15-1-5	2.5/21	Ethylene	5	0.015	straight	Standard
Arterial Pressure Monitoring Set	G58404	APM-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
Arterial Pressure Monitoring Set	G58405	APM-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Set	G58406	APM-3.0-18-1-8	3/20	Ethylene	8	0.018	straight	Standard
Arterial Pressure Monitoring Set	G58407	APM-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
Arterial Pressure Monitoring Set	G58408	APM-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
Arterial Pressure Monitoring Set	G58409	APM-4.0-21-0-6	4/18	Ethylene	6	0.021	straight	Baseplate
Arterial Pressure Monitoring Set	G58410	APM-4.0-21-1-12	4/18	Ethylene	12	0.021	straight	Standard
Arterial Pressure Monitoring Set	G58411	APM-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
Arterial Pressure Monitoring Set - Urethane	G58420	APM-3.0-18-1-8-UR	3/20	Urethane	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Set - Nylon	G58398	APM-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard
Trays								
Arterial Pressure	G58412	APMY-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
Arterial Pressure Monitoring Tray	G58413	APMY-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
Arterial Pressure Monitoring Tray	G58414	APMY-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
Arterial Pressure Monitoring Tray	G58415	APMY-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
Arterial Pressure Monitoring Tray	G58416	APMY-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Tray	G58417	APMY-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
Arterial Pressure Monitoring Tray	G58418	APMY-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
Arterial Pressure Monitoring Tray	G58419	APMY-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
Arterial Pressure Monitoring Tray - Nylon	G58399	APMY-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard



SPECTRUM VASCULAR **Arterial**
PRESSURE MONITORING SETS & TRAYS



For pricing, product specs and additional information on arterial sets & trays, please contact your local Cook Medical or Spectrum Vascular representative.



www.spectrumvascular.com