


# Arterial Pressure Monitoring Sets & Trays

# Arterial Pressure Monitoring Sets & Trays



For pricing, product specs and additional information on arterial sets & trays, please contact your Spectrum Vascular representative.

## We've Got You Covered.

With over 40 years\* as a recognized leader in minimally invasive technologies, you can rely on us. Our portfolio is diverse and is designed to meet the challenges you face.



 **SPECTRUM VASCULAR™**  
spectrumvascular.com

 **SPECTRUM VASCULAR™**  
spectrumvascular.com

Spectrum Vascular and the Spectrum Vascular logo are trademarks of Spectrum Vascular. Spectrum Vascular acquired the vascular access business of Cook Medical (2023) and Angiodynamics (2024). Product continues to be manufactured and labeled / branded as Cook and Angiodynamics pending transition.

© 2026 Spectrum Vascular. All rights reserved. MKT-001 REV. 1

\* Spectrum Vascular acquired the Cook Medical Arterial Sets and Trays portfolio in 2023.

# Arterial Pressure Monitoring Sets & Trays

## Why Spectrum Vascular

Used for arterial pressure monitoring and blood sampling.<sup>†,††,‡</sup>

### FEATURES & BENEFITS

- Inserted via the minimally invasive Seldinger insertion technique to minimize vessel trauma
- Lines available for radial or femoral insertion to increase placement options
- Femoral access set is available with extra-firm nylon catheter

### PRIMARY USE AREAS

- Intensive Care Unit
- Perioperative
- Emergency Room

Available in a range of sizes and design configurations.

Product Name	Order Number	Reference Part Number (RPN)	Catheter Fr/ Equivalent (gauge)	Material	Catheter Length (cm)	Wire Guide Diameter (inch)	Wire Guide Tip Configuration	Needle Baseplate
<b>SETS</b>								
Arterial Pressure Monitoring Set	G58400	APM-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
Arterial Pressure Monitoring Set	G58401	APM-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
Arterial Pressure Monitoring Set	G58402	APM-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
Arterial Pressure Monitoring Set	G58403	APM-2.5-15-1-5	2.5/21	Ethylene	5	0.015	straight	Standard
Arterial Pressure Monitoring Set	G58404	APM-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
Arterial Pressure Monitoring Set	G58405	APM-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Set	G58406	APM-3.0-18-1-8	3/20	Ethylene	8	0.018	straight	Standard
Arterial Pressure Monitoring Set	G58407	APM-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
Arterial Pressure Monitoring Set	G58408	APM-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
Arterial Pressure Monitoring Set	G58409	APM-4.0-21-0-6	4/18	Ethylene	6	0.021	straight	Baseplate
Arterial Pressure Monitoring Set	G58410	APM-4.0-21-1-12	4/18	Ethylene	12	0.021	straight	Standard
Arterial Pressure Monitoring Set	G58411	APM-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
Arterial Pressure Monitoring Set - Urethane	G58420	APM-3.0-18-1-8-UR	3/20	Urethane	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Set - Nylon	G58398	APM-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard
<b>TRAYS</b>								
Arterial Pressure Monitoring Tray	G58412	APMY-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
Arterial Pressure Monitoring Tray	G58413	APMY-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
Arterial Pressure Monitoring Tray	G58414	APMY-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
Arterial Pressure Monitoring Tray	G58415	APMY-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
Arterial Pressure Monitoring Tray	G58416	APMY-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
Arterial Pressure Monitoring Tray	G58417	APMY-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
Arterial Pressure Monitoring Tray	G58418	APMY-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
Arterial Pressure Monitoring Tray	G58419	APMY-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
Arterial Pressure Monitoring Tray - Nylon	G58399	APMY-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard

**~1 in 6**  
individuals hospitalised  
with influenza had  
severe outcomes<sup>1</sup>

Intended for arterial blood pressure monitoring and blood sampling in adult and pediatric patients.

† 2.5 French catheters are intended for patients from birth and older.

†† 3.0 and 4.0 French catheters are intended for patients 1 year and older.

‡ 5.0 French catheters are intended for patients aged 12 years and older.

1. Sumner, et al. Severity of influenza-associated hospitalisations by influenza virus type and subtype in the USA, 2010-19: a repeated cross-sectional study, The Lancet Microbe, Volume 4, Issue 11, 2023, Pages e903-e912, ISSN 2666-5247, [https://doi.org/10.1016/S2666-5247\(23\)00187-8](https://doi.org/10.1016/S2666-5247(23)00187-8).