



# Elevate your arterial line placement.

---

ARTERIAL PRESSURE  
MONITORING CATHETERS





# Pressure Monitoring Catheters

Used for arterial pressure monitoring and blood sampling.

## FEATURES + BENEFITS

- Inserted via the minimally invasive Seldinger insertion technique to minimize vessel trauma.
- Lines available for radial or femoral insertion to increase treatment options.
- Femoral access set available with nylon catheter.

UPN	Reference Part Number	Catheter Fr/ Equivalent gage	Material	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Tip Configuration	Needle Baseplate
<b>Sets</b>							
G58400	APM-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
G58401	APM-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
G58402	APM-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
G58403	APM-2.5-15-1-5	2.5/21	Ethylene	5	0.015	straight	Standard
G58404	APM-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
G58405	APM-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
G58406	APM-3.0-18-1-8	3/20	Ethylene	8	0.018	straight	Standard
G58407	APM-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
G58408	APM-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
G58409	APM-4.0-21-0-6	4/18	Ethylene	6	0.021	straight	Baseplate
G58410	APM-4.0-21-1-12	4/18	Ethylene	12	0.021	straight	Standard
G58411	APM-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
G58420	APM-3.0-18-1-8-UR	3/20	Urethane	5	0.018	straight	Baseplate
G58398	APM-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard
<b>Trays</b>							
G58412	APMY-2.5-15-0-2.5	2.5/21	Ethylene	2.5	0.015	straight	Baseplate
G58413	APMY-2.5-15-0-5	2.5/21	Ethylene	5	0.015	straight	Baseplate
G58414	APMY-2.5-15-1-2.5	2.5/21	Ethylene	2.5	0.015	straight	Standard
G58415	APMY-3.0-18-0-8	3/20	Ethylene	8	0.018	straight	Baseplate
G58416	APMY-3.0-18-0-5	3/20	Ethylene	5	0.018	straight	Baseplate
G58417	APMY-3.0-18-1-5	3/20	Ethylene	5	0.018	straight	Standard
G58418	APMY-4.0-21-0-12	4/18	Ethylene	12	0.021	straight	Baseplate
G58419	APMY-4.0-21J-1-12	4/18	Ethylene	12	0.021	straight & 2 mm J	Standard
G58399	APMY-5.0-35J-1-15-NY	5/16	Nylon	15	0.035	straight & 3 mm J	Standard

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details. Sets include arterial catheter, introducer needle, wire guide, dilator (femoral artery access sets only) and syringe.