

Infusion for emergent situations.

INTRAOSSEROUS
INFUSION NEEDLES





Intraosseous Infusion Needle

Used primarily during pediatric emergencies as an alternative to unsuccessful intravenous access to allow for effective infusion of resuscitative drugs or fluids.

FEATURES + BENEFITS

- Two laterally opposed sideports positioned near the needle cannula's distal tip ensure proper flow if the needle tip is obstructed by posterior bone cortex.
- The positioning mark that is 1 cm proximal to the needle cannula's distal tip provides a visual reference point for depth determination.

UPN	Reference Part Number	Needle Gage	Cannula / Needle Length cm
G05312	C-DIN-14-3.0-T45-DKM	14	3
G05986	C-DIN-16-2.5-T45-DKM	16	2.5
G04958	C-DIN-16-3.0-T45-DKM	16	3
G04954	C-DIN-16-4.0-T45-DKM	16	4
G07004	C-DIN-18-2.5-T45-DKM	18	2.5
G05133	C-DIN-18-3.0-T45-DKM	18	3

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details.