

Making a material difference.

PERIPHERALLY INSERTED
CENTRAL CATHETERS (PICC)






BioFlo® with Endexo® Technology

The BioFlo® PICC Catheter is the first - and only - PICC catheter with Endexo® Technology, providing a catheter material more resistant to in-vitro thrombus accumulation, compared to commonly used catheters (based on platelet count).¹

Endexo® Technology is a permanent and non-eluting polymer that is “blended” into the polyurethane from which the catheter is made. It is present throughout the catheter, including the intraluminal, extraluminal, and cut surface of the catheter. Endexo® Technology remains present for the life of the catheter.

● PICCs

Peripherally inserted central catheters are indicated for short and long-term peripheral access to the central venous system for intravenous therapy.



Learn more about BioFlo®
technology, including links
to research on BioFlo®, at
[spectrumvascular.com](https://www.spectrumvascular.com)

1. The reduction in thrombus accumulation (based on platelet count) is supported by acute in-vitro testing. Pre-clinical in-vitro evaluations do not necessarily predict clinical performance with respect to thrombus formation.

BioFlo® PICC CATHETER ONLY KITS

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
92cm Tape Measure
End Cap(s)
Catheter Securement Device

BioFlo® PICC MST KITS WITH 70CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
70 cm Hydrophilic-Coated Guidewire
Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

BioFlo® PICC MST KITS WITH 70CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
70 cm Hydrophilic-Coated Guidewire
Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

BioFlo® PICC INTERMEDIATE KITS

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Face Mask
2 - Tourniquets
2 - 92 cm Tape Measures
Scissors
Absorbant Poly-Lined Drape (24" x 36")
Absorbant Poly-Lined Drape (24" x 38")
3 mL Chloraprep
3 mL Luer Lock Syringe
Safety Hypodermic Needle (25G x 0.625")
45 cm Guidewire
Safety Scalpel
Sheath/Dilator Introducer
Non-Serated Forceps
Gauze Pads (2" x 2")
Gauze Pads (4" x 4")
10 mL Luer Lock Syringe(s)
Saline Ampule(s)
Safety Ampule Breaker
Filter Straw
Saline Label(s)
End Cap(s)
Catheter Securement Device
Transparent Dressing
CSR Wrap

^Lidocaine based on kit configuration

BioFlo® PICC MAXIMAL BARRIER KITS

Catheter
Gown (XL)
Saline - Pre-filled 10 mL Syringes
Transducer Probe Cover Set
Maximal Barrier Body Drape(s)*
Underarm Drape
Port Protector
Tuohy Borst with Side Port
Stiffening Stylet or 70 cm Nitinol Guidewire**
Gauze Pads (4" x 4")
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
Sharps container
Guidewire 0.018" x 45 cm Nitinol
Safety Scalpel
Sheath/Dilator Introducer
3 mL Clear Chloraprep
3 mL Orange Tinted Chloraprep
Statlock Stabilization Device
Transparent Dressing
Absorbant Towel (14" x 25")
Tourniquet
Surgical Tape
CSR Wrap
Mask with Earloops
Hair Net
Gloves
Safety Ampule Breaker
Lidocaine HCL 1%
5 mL Ampule
Lidocaine Label
5 mL Syringe
Filter Straw
Needle (25G x 0.625")
Tape Measure
Scissors
Patient Checklist
Stop Sign

*2 options:
1) Full Body Drape (70" x 100") and Fenestrated Drape
2) Dual Fenestrated Drape (73" x 114")

**Based on kit configuration

BioFlo® PICCs • Ordering Guide

BIOFLO® POWER INJECTABLE PICC – NON-VALVED

STAINLESS STEEL GUIDEWIRE

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
-----	-----------------------	--------	--------------------	----------------------------	----------------------

CATHETER KITS

H965459540	45-954	Single	3	4	20
H965458750	45-875	Single	4	6	17
H965458700	45-870	Dual	5	7	17.5/17.5
H965459620	45-962	Triple	6	7	16.5/19/19

MST® KITS WITH 70 CM WIRE

H965459550	45-955	Single	3	4	20
H965458760	45-876	Single	4	6	17
H965458570	45-857	Single	5	6	15.5
H965458710	45-871	Dual	5	7	17.5/17.5
H965458670	45-867	Dual	6	7	16.5/16.5
H965459630	45-963	Triple	6	7	16.5/19/19

IR KITS WITH 145 CM WIRE

H965458770	45-877	Single	4	6	17
H965458580	45-858	Single	5	6	15.5
H965458720	45-872	Dual	5	7	17.5/17.5
H965458680	45-868	Dual	6	7	16.5/16.5
H965459640	45-964	Triple	6	7	16.5/19/19

INTERMEDIATE MST KITS WITH 45 CM WIRE

H965458730	45-873	Dual	5	7	17.5/17.5
H965458690	45-869	Dual	6	7	16.5/16.5

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details.

BioFlo® PICCs • Ordering Guide

BIOFLO® POWER INJECTABLE PICC – NON-VALVED (CONTINUED)

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
INTERMEDIATE MST KITS WITH 45 CM WIRE AND LIDOCAINE					
H965459580	45-958	Single	3	4	20
H965459660	45-966	Triple	6	7	16.5/19/19
NITINOL GUIDEWIRE					
MST KITS WITH 70 CM WIRE					
H965459590	45-959	Single	3	4	20
H965458800	45-880	Single	4	6	17
H965458960	45-896	Single	5	6	15.5
H965458950	45-895	Dual	5	7	17.5/17.5
H965459680	45-968	Triple	6	7	16.5/19/19
IR KITS WITH 145 CM WIRE					
H965458850	45-885	Single	4	6	17
H965458860	45-886	Single	5	6	15.5
H965458870	45-887	Dual	5	7	17.5/17.5
H965458880	45-888	Dual	6	7	16.5/16.5
H965459670	45-967	Triple	6	7	16.5/19/19

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details.

BioFlo® PICCs • Ordering Guide

MAXIMAL BARRIER NURSING KIT

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
MAXIMAL BARRIER NURSING KITS PG					
H965750231	75-023	Single	3	4	20
H965750121	75-012	Single	4	6	17
H965750131	75-013	Single	5	6	15.5
H965750141	75-014	Dual	5	7	17.5/17.5
H965750151	75-015	Dual	6	7	16.5/16.5
H965750241	75-024	Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH 70 CM NITINOL WIRE

H965750315	75-031	Single	4	6	17
H965750325	75-032	Single	5	6	15.5
H965750335	75-033	Dual	5	7	17.5/17.5
H965750345	75-034	Dual	6	7	16.5/16.5
H965750355	75-035	Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH DUAL FENESTRATED DRAPE

H965750541	75-054	Single	3	4	20
H965750671	75-067	Single	4	6	17
H965750441	75-044	Single	5	6	15.5
H965750451	75-045	Dual	5	7	17.5/17.5
H965750461	75-046	Dual	6	7	16.5/16.5
H965750551	75-055	Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH 70 CM NITINOL WIRE AND DUAL FENESTRATED DRAPE

H965750561	75-056	Single	4	6	17
H965750571	75-057	Single	5	6	15.5
H965750581	75-058	Dual	5	7	17.5/17.5
H965750591	75-059	Dual	6	7	16.5/16.5
H965750601	75-060	Triple	6	7	16.5/16.5

BioFlo® PICCs • Ordering Guide

BIOFLO® POWER INJECTABLE PICC – W/ PASV VALVE TECHNOLOGY

STAINLESS STEEL GUIDEWIRE

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
CATHETER KITS					
H965458010	45-801	Single	3	4	20
H965458110	45-811	Single	4	6	17
H965458310	45-831	Dual	5	7	17.5/17.5
H965458410	45-841	Triple	6	7	16.5/19/19
INTERMEDIATE MST KITS WITH 45 CM WIRE					
H965458040	45-804	Single	3	4	20
INTERMEDIATE MST KITS WITH 45 CM WIRE AND LIDOCAINE					
H965458050	45-805	Single	3	4	20
H965458150	45-815	Single	4	6	17
H965458350	45-835	Dual	5	7	17.5/17.5
H965458450	45-845	Triple	6	7	16.5/19/19
MST KITS WITH 70 CM WIRE					
H965458020	45-802	Single	3	4	20
H965458120	45-812	Single	4	6	17
H965458170	45-817	Single	5	6	15.5
H965458320	45-832	Dual	5	7	17.5/17.5
H965458370	45-837	Dual	6	7	16.5/16.5
H965458420	45-842	Triple	6	7	16.5/19/19
IR KITS WITH 145 CM WIRE					
H965458130	45-813	Single	4	6	17
H965458180	45-818	Single	5	6	15.5
H965458330	45-833	Dual	5	7	17.5/17.5
H965458380	45-838	Dual	6	7	16.5/16.5
H965458430	45-843	Triple	6	7	16.5/19/19

BioFlo® PICCs • Ordering Guide

BIOFLO® POWER INJECTABLE PICC – W/ PASV VALVE TECHNOLOGY (CONTINUED)

NITINOL GUIDEWIRE

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
MST KITS WITH 70 CM WIRE					
H965458940	45-894	Single	3	4	20
H965458910	45-891	Single	4	6	17
H965458920	45-892	Single	5	6	15.5
H965458890	45-889	Dual	5	7	17.5/17.5
H965458900	45-890	Dual	6	7	16.5/16.5
H965458930	45-893	Triple	6	7	16.5/19/19

IR KITS WITH 145 CM WIRE

H965458810	45-881	Single	4	6	17
H965458820	45-882	Single	5	6	15.5
H965458830	45-883	Dual	5	7	17.5/17.5
H965458840	45-884	Triple	6	7	16.5/19/19

HYBRID PICC KITS

CATHETER KIT

H965459410	45-941	Hybrid Triple	6	7	16.5/19/19
------------	--------	---------------	---	---	------------

MST KIT WITH 70 CM WIRE

H965459420	45-942	Hybrid Triple	6	7	16.5/19/19
------------	--------	---------------	---	---	------------

MST KIT WITH 70 CM NITINOL WIRE

H965458970	45-897	Hybrid Triple	6	7	16.5/19/19
------------	--------	---------------	---	---	------------

IR KIT WITH 145 CM WIRE

H965459430	45-943	Hybrid Triple	6	7	16.5/19/19
------------	--------	---------------	---	---	------------

BioFlo® PICCs • Ordering Guide

BIOFLO® POWER INJECTABLE PICC – W/ PASV VALVE TECHNOLOGY (CONTINUED)

MAXIMAL BARRIER NURSING KIT

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
MAXIMAL BARRIER NURSING KITS PG					
H965750161	75-016	Single	3	4	20
H965750171	75-017	Single	4	6	17
H965750181	75-018	Single	5	6	15.5
H965750191	75-019	Dual	5	7	17.5/17.5
H965750211	75-021	Triple	6	7	16.5/19/19
H965750221	75-022	Hybrid Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH 70 CM NITINOL WIRE

H965750365	75-036	Single	4	6	17
H965750375	75-037	Single	5	6	15.5
H965750385	75-038	Dual	5	7	17.5/17.5
H965750415	75-041	Triple	6	7	16.5/19/19
H965750425	75-042	Hybrid Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH DUAL FENESTRATED DRAPE

H965750471	75-047	Single	3	4	20
H965750481	75-048	Single	4	6	17
H965750491	75-049	Single	5	6	15.5
H965750501	75-050	Dual	5	7	17.5/17.5
H965750511	75-051	Dual	6	7	16.5/16.5
H965750521	75-052	Triple	6	7	16.5/19/19
H965750531	75-053	Hybrid Triple	6	7	16.5/19/19

MAXIMAL BARRIER NURSING KITS WITH 70 CM NITINOL WIRE AND DUAL FENESTRATED DRAPE

H965750611	75-061	Single	4	6	17
H965750621	75-062	Single	5	6	15.5
H965750631	75-063	Dual	5	7	17.5/17.5
H965750641	75-064	Dual	6	7	16.5/16.5
H965750651	75-065	Triple	6	7	16.5/19/19
H965750661	75-066	Hybrid Triple	6	7	16.5/19/19

Xcela Power Injectable PICC CATHETER KITS

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
92cm Tape Measure
End Cap(s)
Catheter Securement Device

Xcela Power Injectable PICC MST KITS WITH 70CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
70 cm Hydrophilic-Coated Guidewire

Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

Xcela Power Injectable PICC IR KITS WITH 145CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
145 cm Hydrophilic-Coated Guidewire

Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

Xcela Power Injectable PICC MAXIMAL BARRIER KITS

Catheter
Gown (XL)
Saline - Pre-filled 10 mL Syringes
Transducer Probe Cover Set
Maximal Barrier Body Drape(s)*
Underarm Drape
Port Protector
Tuohy Borst with Side Port
Stiffening Stylet or 70 cm Nitinol Guidewire**
Gauze Pads (4" x 4")
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
Sharps container
Guidewire 0.018" x 45 cm Nitinol
Safety Scalpel
Sheath/Dilator Introducer
3 mL Clear Chloraprep
3 mL Orange Tinted Chloraprep
Statlock Stabilization Device
Transparent Dressing
Absorbant Towel (14" x 25")
Tourniquet
Surgical Tape
CSR Wrap
Mask with Earloops
Hair Net
Gloves
Safety Ampule Breaker
Lidocaine HCL 1%
5 mL Ampule
Lidocaine Label
5 mL Syringe
Filter Straw
Needle (25G x 0.625")
Tape Measure
Scissors
Patient Checklist
Stop Sign

*2 options:
1) Full Body Drape (70" x 100") and Fenestrated Drape
2) Dual Fenestrated Drape (73" x 114")

**Based on kit configuration

Xcela Power Injectable PICCs • Ordering Guide

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)	Guidewire
CATHETER KITS						
M001456750†	45-675	Single	4	6	17	Stainless Steel
M001456700†	45-670	Dual	5	7	17.5/17.5	Stainless Steel
SAFETY MST KITS WITH 60 CM WIRE						
M001456520	45-652	Single	4	6	17	Stainless Steel
M001456760	45-676	Single	4	6	17	Stainless Steel
M001456570	45-657	Single	5	6	15.5	Stainless Steel
M001456710	45-671	Dual	5	7	17.5/17.5	Stainless Steel
M001456670	45-667	Dual	6	7	16.5/16.5	Stainless Steel
SAFETY IR KITS WITH 145 CM WIRE						
M001456770	45-677	Single	4	6	17	Stainless Steel
M001456580	45-658	Single	5	6	15.5	Stainless Steel
M001456720	45-672	Dual	5	7	17.5/17.5	Stainless Steel
M001456680	45-668	Dual	6	7	16.5/16.5	Stainless Steel
INTERMEDIATE SAFETY MST KITS WITH 45 CM WIRE						
M001456780	45-678	Single	4	6	17	Stainless Steel
M001456640†	45-664	Dual	5	7	17.5/17.5	Stainless Steel
M001456730	45-673	Dual	5	7	17.5/17.5	Stainless Steel

Xcela Power Injectable PICCs • Ordering Guide

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)	Guidewire
MST KITS WITH 70 CM WIRE						
H965251010	25-101	Single	4	6	17	Nitinol
H965251030	25-103	Single	5	6	15.5	Nitinol
H965251050	25-105	Dual	5	7	17.5/17.5	Nitinol
H965251080	25-108	Dual	6	7	16.5/16.5	Nitinol

IR KITS WITH 145 CM WIRE						
H965251020	25-102	Single	4	6	17	Nitinol
H965251040	25-104	Single	5	6	15.5	Nitinol
H965251060	25-106	Dual	5	7	17.5/17.5	Nitinol
H965251090	25-109	Dual	6	7	16.5/16.5	Nitinol

MAXIMAL BARRIER NURSING KITS						
H965750011	75-001	Single	4	6	17	Nitinol
H965750021	75-002	Single	5	6	15.5	Nitinol
H965750031	75-003	Dual	5	7	17.5/17.5	Nitinol
H965750041 ¹	75-004	Dual	6	7	16.5/16.5	Nitinol

MAXIMAL BARRIER WITH DUAL FENESTRATED DRAPE KITS						
H965750681	75-068	Single	4	6	17	Nitinol
H965750691	75-069	Single	5	6	15.5	Nitinol
H965750701	75-070	Dual	5	7	17.5/17.5	Nitinol
H965750711	75-071	Dual	6	7	16.5/16.5	Nitinol

UPN	Reference Part Number	Proximal/Distal Outer Diameter	Total Length	Guidewire/Tip Material	Units
SINGLE STERILE NITINOL HYDROPHILIC-COATED GUIDEWIRE					
H965457461	45-746	0.018"	70 cm	Nitinol/Tungsten	Box/3
H965457471	45-747	0.018"	145 cm	Nitinol/Tungsten	Box/3

Xcela PICC with PASV[®] Valve Technology CATHETER KITS

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
92cm Tape Measure
End Cap(s)
Catheter Securement Device

Xcela PICC with PASV[®] Valve Technology MST KITS WITH 70CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
70 cm Hydrophilic-Coated Guidewire

Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

Xcela PICC with PASV[®] Valve Technology IR KITS WITH 145CM WIRE

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Safety Introducer Needle (21G x 2.75")
Standard Introducer Needle (21G x 2.75")
145 cm Hydrophilic-Coated Guidewire

Safety Scalpel
10 cm Sheath/Dilator
92 cm Tape Measure
10 mL Luer Lock Syringe(s)
End Cap(s)
Catheter Securement Device

Xcela Hybrid PICC with PASV[®] Valve Technology INTERMEDIATE KITS

Catheter
Wire Guide/Flush Assembly
PTFE Coated Siffening Wire
Face Mask
2 - Tourniquets
2 - 92 cm Tape Measures
Scissors
Absorbant Poly-Lined Drape (24" x 36")
Absorbant Poly-Lined Drape (24" x 38")
3 mL Chloraprep
3 mL Luer Lock Syringe
Safety Hypodermic Needle (25G x 0.625")
45 cm Guidewire
Safety Scalpel
Sheath/Dilator Introducer
Non-Serated Forceps
Gauze Pads (2" x 2")
Gauze Pads (4" x 4")
10 mL Luer Lock Syringe(s)
Saline Ampule(s)
Safety Ampule Breaker
Filter Straw
Saline Label(s)
End Cap(s)
Catheter Securement Device
Transparent Dressing
CSR Wrap

^Lidocaine based on kit configuration

Xcela PICC with PASV* Valve Technology • Ordering Guide

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)	Guidewire
CATHETER KITS						
H965457110	45-711	Single	4	6	17	Stainless Steel
H965457160	45-716	Single	5	6	15.5	Stainless Steel
H965457310	45-731	Dual	5	7	17.5/17.5	Stainless Steel
H965457410†	45-741	Triple	6	7	16.5/19/19	Stainless Steel
INTERMEDIATE MST KITS WITH 45 CM WIRE						
H965457140	45-714	Single	4	6	17	Stainless Steel
H965457340	45-734	Dual	5	7	17.5/17.5	Stainless Steel
H965457440	45-744	Triple	6	7	16.5/19/19	Stainless Steel
INTERMEDIATE MST KITS WITH 45 CM WIRE AND LIDOCAINE						
H965457150	45-715	Single	4	6	17	Stainless Steel
H965457350	45-735	Dual	5	7	17.5/17.5	Stainless Steel
H965457450	45-745	Triple	6	7	16.5/19/19	Stainless Steel
MST KITS WITH 70 CM WIRE						
H965457120	45-712	Single	4	6	17	Stainless Steel
H965457170	45-717	Single	5	6	15.5	Stainless Steel
H965457320	45-732	Dual	5	7	17.5/17.5	Stainless Steel
H965457370†	45-737	Dual	6	7	16.5/16.5	Stainless Steel
H965457420	45-742	Triple	6	7	16.5/19/19	Stainless Steel
IR KITS WITH 145 CM WIRE						
H965457130	45-713	Single	4	6	17	Stainless Steel
H965457180	45-718	Single	5	6	15.5	Stainless Steel
H965457330	45-733	Dual	5	7	17.5/17.5.5	Stainless Steel
H965457380	45-738	Dual	6	7	16.5/16.5	Stainless Steel
H965457430	45-743	Triple	6	7	16.5/19/19	Stainless Steel

Xcela PICC with PASV* Valve Technology • Ordering Guide

UPN	Reference Part Number	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)	Guidewire
MST KITS WITH 70 CM WIRE						
H965251210	25-121	Single	4	6	17	Nitinol
H965251230	25-123	Single	5	6	15.5	Nitinol
H965251250	25-125	Dual	5	7	17.5/17.5	Nitinol
H965251290	25-129	Triple	6	7	16.5/19/19	Nitinol
IR KITS WITH 145 CM WIRE						
H965251220	25-122	Single	4	6	17	Nitinol
H965251240	25-124	Single	5	6	15.5	Nitinol
H965251260	25-126	Dual	5	7	17.5/17.5	Nitinol
H965251300†	25-130	Triple	6	7	16.5/19/19	Nitinol

Xcela Hybrid PICC with PASV Valve Technology • Ordering Guide

UPN	Reference Part Number	Description	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
INTERMEDIATE MST KITS						
H965952440	95-244	Intermediate MST Kits with 45 cm Wire	3	6	7	16.5/19/19
H965952450	95-245	Intermediate MST Kits with 45 cm Wire and Lidocaine	3	6	7	16.5/19/19
MST KITS						
H965952420	95-242	MST [†] Kits with 70 cm Wire	3	6	7	16.5/19/19
H965952460 [†]	95-246	MST [†] Kits with 70 cm Nitinol Wire	3	6	7	16.5/19/19
IR KITS						
H965952430	95-243	IR Kits with 145 cm Wire	3	6	7	16.5/19/19
H965952470	95-247	IR Kits with 145 cm Nitinol Wire	3	6	7	16.5/19/19

MST = Modified Seldinger Technique


[†] These items are available Made to Order. Contact Spectrum Vascular Customer Service or your local Territory Manager for details concerning these products.

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details.

CT, Non-Valved Catheters

The Pro-PICC CT Catheter reduces patient trauma and risks associated with secondary venous access by allowing for the combination of CT Power Injection capabilities and infusion therapy. Flow rates are easily identifiable on the ID tags, which also provide strength and durability for repeated use.

PRO-PICC^{CT}



Printed Identification and Purple Color
denote maximum flow rates for contrast injection

Easy-to-Read Depth Marks
for easy identification of catheter length

Pro-PICC CT MST Kits

PICC (w/ stylet)
Peelable sheath introducer
Mini scalpel
Guidewire
10 cc Syringe
Tape measure(s)

Introducer needle
Needle free valve(s)
Catheter securement device
Patient information pack
Patient chart sticker*
Patient ID card*

Pro-PICC CT IR Kits

PICC (w/ stylet)
Peelable sheath introducer
Mini scalpel
Guidewire
10 cc Syringe
Tape measure(s)

Introducer needle
Needle free valve(s)
Catheter securement device
Patient information pack
Patient chart sticker*
Patient ID card*

Pro-PICC CT Basic Nursing Kits

PICC (w/ stylet)
Peelable sheath introducer
Guidewire
10 cc Syringe
Tape measure(s)
Needle free valve(s)

Catheter securement device
Patient information pack
Patient chart sticker*
Patient ID card*
Safety scalpel
Safety introducer needle

Pro-PICC CT Maximal Barrier Kits

PICC (w/ stylet)
Peelable sheath introducer
Guidewire
Tape measure(s)
Introducer needle
Needle free valve(s)
Catheter securement device
Patient information pack
Patient chart sticker*
Patient ID card*
Safety scalpel
Safety introducer needle
Safety Needle
Tourniquet
5 cc Syringe
Surgical Tape
Pre-filled saline syringe(s)
Chlorhexidine
Tinted chlorhexidine
5 cc Lidocaine
Adhesive wound dressing
Needle stick pad
4" x 4" Gauze
Surgical gloves
Filter straw
Scissors
Steri-strips
Face mask(s)
2" x 2" Gauze
Fenestrated drape
Drape
Safety hypodermic needle
Lidocaine/saline stickers
Full body drape
XL gown
Surgical face mask with eye shield
Bouffant cap
Probe cover with gel

Pro-PICC CT • Ordering Guide


UPN	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
MST KITS				
MD17033101	Single	3	4	19
MD17034101	Single	4	5.5	18
MD17035101	Single	5	7	18
MD17035201	Dual	5	7	18/18
MD17036201	Dual	6	7	18/18
MD17036301	Triple	6	7.5	17/19/19
IR KITS				
MD17033102	Single	3	4	19
MD17034102	Single	4	5.5	18
MD17035102	Single	5	7	18
MD17035202	Dual	5	7	18/18
MD17036202	Dual	6	7	18/18
MD17036302	Triple	6	7.5	17/19/19
BASIC NURSING KITS				
MD17033105	Single	3	4	19
MD17034105	Single	4	5.5	18
MD17035105	Single	5	7	18
MD17035205	Dual	5	7	18/18
MD17036205	Dual	6	7	18/18
MD17036305	Triple	6	7.5	17/19/19
MAXIMAL BARRIER KITS				
MR17033106MB	Single	3	4	19
MR17034106MB	Single	4	5.5	18
MR17035106MB	Single	5	7	18
MR17035206MB	Dual	5	7	18/18
MR17036306MB	Triple	6	7.5	17/19/19

Non-CT, Non-Valved Catheters

Vascu-PICC and Vascu-PICC TL catheters provide a range of options to meet the unique needs of your patients. Vascu-PICC offers a tapered design that may help reduce bleeding at the insertion site, while Vascu-PICC TL features a taperless lumen, ideal for smaller vasculature.



Taperless (Vascu-PICC TL) or
Tapered (Vascu-PICC) Design



Polyurethane Material

Radiopaque
Identifiable
under X-Ray

**Easy-to-Read
Depth Marks**

Vascu-PICC and Vascu-PICC TL MST Kits

PICC (w/ stylet)
 Peelable sheath introducer
 Mini scalpel
 Guidewire
 10 cc Syringe
 Tape measure(s)
 Introducer needle
 Needle free valve(s)
 Catheter securement device
 Patient information pack
 Patient chart sticker*
 Patient ID card*

Vascu-PICC and Vascu-PICC TL IR Kits

PICC (w/ stylet)
 Peelable sheath introducer
 Mini scalpel
 Guidewire
 10 cc Syringe
 Tape measure(s)
 Introducer needle
 Needle free valve(s)
 Catheter securement device
 Patient information pack
 Patient chart sticker*
 Patient ID card*

Vascu-PICC and Vascu-PICC TL Basic Nursing Kits

PICC (w/ stylet)
 Peelable sheath introducer
 Guidewire
 10 cc Syringe
 Tape measure(s)
 Needle free valve(s)
 Catheter securement device
 Patient information pack
 Patient chart sticker*
 Patient ID card*
 Safety scalpel
 Safety introducer needle

Vascu-PICC and Vascu-PICC TL Maximal Barrier Kits

PICC (w/ stylet)
 Peelable sheath introducer
 Guidewire
 Tape measure(s)
 Introducer needle
 Needle free valve(s)
 Catheter securement device
 Patient information pack
 Patient chart sticker*
 Patient ID card*
 Safety scalpel
 Safety introducer needle
 Safety Needle
 Tourniquet
 5 cc Syringe
 Surgical Tape
 Pre-filled saline syringe(s)
 Chlorhexidine
 Tinted chlorhexidine
 5 cc Lidocaine
 Adhesive wound dressing
 Needle stick pad
 4" x 4" Gauze
 Surgical gloves
 Filter straw
 Scissors
 Steri-strips
 Face mask(s)
 2" x 2" Gauze
 Fenestrated drape
 Drape
 Safety hypodermic needle
 Lidocaine/saline stickers
 Full body drape
 XL gown
 Surgical face mask with eye shield
 Bouffant cap
 Probe cover with gel

Vascu-PICC • Ordering Guide

UPN	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
MST KITS				
MD17013101	Single	3	7.5	19
MD17014101	Single	4	4	18
MD17015101	Single	5	6	16
MD17015201	Dual	5	7	18/18
MD17016201*	Dual	6	7	17/17
IR KITS				
MD17015202	Dual	5	7	18/18
BASIC NURSING KITS				
MD17013105	Single	3	7	19
MD17014105	Single	4	4	18
MD17014205	Dual	4	6	19/19
MAXIMAL BARRIER KITS				
MR17014106MB	Single	4	6	18
MR17015206MB	Dual	5	6	18/18

* Spectrum Vascular does not keep high inventory levels of this code. Please check with your sales representative for availability and lead times.


Vascu-PICC TL (Taperless) • Ordering Guide

UPN	Lumens	Outer Diameter (F)	Reverse Taper Diameter (F)	Inner Diameter (Ga.)
BASIC MST KITS				
MD81013101	Single	3	7	19
MD81014101*	Single	4	N/A	18
MD81014201*	Dual	4	N/A	19/19
MD81015201*	Dual	5	N/A	18/18
MAXIMAL BARRIER KITS				
MR81013106MB	Single	3	N/A	19
MR81014106MB	Single	4	N/A	18
MR81014206MB	Dual	4	N/A	19/19
MR81015206MB	Dual	5	N/A	18/18

* Spectrum Vascular does not keep high inventory levels of this code. Please check with your sales representative for availability and lead times.

Vascu-PICC TLS Pediatric

Designed to provide exceptional treatment for your smallest patients. Vascu-PICC TLS features dual and single lumen configurations, a polyurethane body, and high flow rates in a small French size. Since smaller veins require greater precision, it makes sense to offer you more choices in PICC sizes.



Vascu-PICC TLS Pediatric CATHETER ONLY KITS

Vascu-PICC TLS with pre-loaded stylet
Patient information booklet
Needle free valve(s)
Patient chart sticker
Securement device

Vascu-PICC TLS Pediatric BASIC KITS (NO MST)

Vascu-PICC TLS with pre-loaded stylet
Needle free valve(s)
Securement device
10 cc Syringe
2" x 2" Gauze
Tape measure(s)
Tourniquet

Vascu-PICC TLS Pediatric MAXIMAL BARRIER KITS

2" x 2" Gauze
Tape measure(s)
Safety needle
Pre-filled saline syringes (3)
Safety scalpel
Scissors
Chlorhexidine (2)
Surgical tape
37" x 24" Fenestrated drape
Steri-strips
Adhesive dressing
Needle stick pad
5 cc Lidocaine HCL 1%
Filter straw
Lidocaine/saline stickers
Surgical gloves (3)
40" x 70" Drape
Face masks (2)
36" x 24" Drape
5 cc Syringe
Bouffant cap
4" x 4" Gauze
XL Gown

Vascu-PICC TLS Pediatric BASIC OTN KITS

Vascu-PICC TLS with pre-loaded stylet
Needle free valve(s)
Patient chart sticker
Securement device
10 cc Syringe
2" x 2" Gauze
Tape measure(s)
Tearaway introducer
Tourniquet

Vascu-PICC TLS Pediatric BASIC MST KITS

Vascu-PICC TLS with pre-loaded stylet
Needle free valve(s)
Securement device
10 cc Syringe
2" x 2" Gauze
Tape measure(s)
Tearaway introducer
Tourniquet
Guidewire
Safety needle

Vascu-PICC TLS Pediatric • Ordering Guide

UPN	Lumens	Outer Diameter (F)	Length (cm)	Inner Diameter (Ga.)
CATHETER ONLY KITS				
MD17011100	Single	1.9	50	23
MD17012600	Dual	2.6	50	21/23
BASIC OTN KITS				
MD17011121	Single	1.9	20	23
MD17011101	Single	1.9	50	23
MD17012621	Dual	2.6	20	21/23
MD17012601	Dual	2.6	50	21/23
BASIC KITS (NO INSERTION COMPONENTS)				
MD17021121	Single	1.9	20	23
MD17022601	Dual	2.6	50	21/23
BASIC MST KITS				
MD170111224S	Single	1.9	20	23
MD170111024S	Single	1.9	50	23
MD170126224S	Dual	2.6	20	21/23
MD170126024S	Dual	2.6	50	21/23
MAXIMAL BARRIER KITS				
MRVSP19001	Single	1.9	20	23
MRVSP26001	Dual	2.6	50	21/23
ACCESSORIES				
MR17019001	Pediatric PICC Procedure w/ Maximal Barrier			

Some products or part numbers may not be available in all markets. Contact your local Spectrum Vascular sales representative or Customer Service for details.